

## **SAMPLE PREPARATION GUIDE FOR FRESH FROZEN SAMPLES ON Stereo-seq CHIPS OR Stereo-seq CHIP SLIDES**

### USER MANUAL



# REVISION HISTORY

**Manual Version:** A  
**Date:** Oct. 2024  
**Description:** Initial release.

---

**Manual Version:** B  
**Date:** Aug. 2025  
**Description:**

- Removed Chapter 3: CRYOSECTION PREPARATION. The instructions for cryosection preparation and tissue mounting have been integrated into the respective user manuals of each kit.
- Added information and sample preparation requirements for the Stereo-seq Transcriptomics Large Chip Designs Solution.

**Note: Please download the latest version of the manual and use it with the corresponding Stereo-seq Permeabilization kit.**

©2025 STOmics Tech Co., Ltd.

All rights reserved.

1. The products shall be for research use only, not for use in diagnostic procedures.
2. The Content on this manual may be protected in whole or in part by applicable intellectual property laws. STOmics Tech Co., Ltd. and/or corresponding right subjects own their intellectual property rights according to law, including but not limited to trademark rights, copyrights, etc.
3. STOmics Tech Co., Ltd. does not grant or imply the right or license to use any copyrighted content or trademark (registered or unregistered) of ours or any third party's. Without our written consent, no one shall use, modify, copy, publicly disseminate, change, distribute, or publish the program or contents of this manual without authorization, and shall not use the design or the design skills to use or take possession of the trademarks, the logo, or other proprietary information (including images, text, web design or form) of ours or those of our affiliates.
4. Nothing contained herein is intended to or shall be construed as any warranty, expression or implication of the performance of any products listed or described herein. Any and all warranties applicable to any products listed herein are set forth in the applicable terms and conditions of sale accompanying the purchase of such product. STOmics Tech Co., Ltd. makes no warranty and hereby disclaims any and all warranties as to the use of any third-party products or protocols described herein.

# TABLE OF CONTENTS



## CHAPTER 1: INTRODUCTION

1.1. Introduction	1
-------------------	---

## CHAPTER 2: SAMPLE PREPARATION

2.1. Sample Requirements for Fresh Frozen Tissue	4
2.2. Sample Embedding for Fresh Frozen Tissue	7
2.3. Sample Storage and Transportation	11



**NOTE:** Additional operation tips and guidance.




**QUALITY CHECK POINT**



**CAUTION:** Proceed with extra care; improper handling or carelessness may cause experimental failure or accidents.

# CHAPTER 1

## INTRODUCTION



## 1.1. Introduction

The STOmics Stereo-seq Transcriptomics Set is intended for generating a spatially-resolved 3' mRNA library from biological tissue sections requires a Stereo-seq Chip Slide with intact tissue sections as input. Prior to the Stereo-seq Transcriptomics workflow, optimal permeabilization time needs to be determined using the Stereo-seq Permeabilization Set for fresh frozen tissues. This guide provides general guidelines on how to properly perform tissue embedding, sectioning and mounting to better preserve the morphological quality of the tissue sections and the integrity of mRNA transcripts.

The Stereo-seq Chip Slides or Stereo-seq Chips prepared in this guideline are part of these products:

Stereo-seq Permeabilization Set for Chip-on-a-slide V1.1,  
Cat. No.: 211SP11118 / 211SP11118-CG

Stereo-seq Transcriptomics Set for Chip-on-a-slide V1.3,  
Cat. No.: 211ST13114 / 211ST13114-CG

Stereo-seq Transcriptomics Set for Chip-on-a-slide V1.3 (0.5cm \* 0.5cm),  
Cat. No.: 211ST13004 / 211ST13004-CG



Stereo-seq Permeabilization Set (1cm \* 2cm) V1.1, Cat. No.: 111SP11124 / 111SP11124-CG

Stereo-seq Transcriptomics Set (1cm \* 2cm) V1.3, Cat. No.: 111ST13122 / 111ST13122-CG

Stereo-seq Permeabilization Set (2cm \* 2cm) V1.1, Cat. No.: 111SP11224 / 111SP11224-CG

Stereo-seq Transcriptomics Set (2cm \* 2cm) V1.3, Cat. No.: 111ST13221 / 111ST13221-CG


Stereo-seq Permeabilization Set (2cm \* 3cm) V1.1, Cat. No.: 111SP11234 / 111SP11234-CG

Stereo-seq Transcriptomics Set (2cm \* 3cm) V1.3, Cat. No.: 111ST13231 / 111ST13231-CG



### Stereo-seq Chip P Slide and Chip T Slide

- Includes two Stereo-seq Chip P Slides containing **four** Chip P (1cm \* 1cm) on each slide in each Stereo-seq Chip P Kit.
- Includes four Stereo-seq Chip T Slides containing **one** Chip T (1cm \* 1cm) on each slide or **one** Chip T (0.5cm \* 0.5cm) on each slide in each Stereo-seq Chip T Kit.
- Stereo-seq Chip P Slides and Stereo-seq Chip T Slides are differentiated by a laser-engraved label at the end of the slide.



### Stereo-seq Chip P and Chip T (for Large Chip Designs Solution)

- Includes **four** Stereo-seq Chip P in each Stereo-seq Chip P Kit.
- Includes **two** Stereo-seq Chip T (1cm \* 2cm) or **one** Stereo-seq Chip T (2cm \* 2cm) or **one** Chip T (2cm \* 3cm) in in each Stereo-seq Chip T Kit.

### Stereo-seq Chip (Slide) Storage

Always store unused slides in their original slide container and keep them sealed with a sealable aluminum bag at 2°C ~ 8°C. Keep sealed with tape or a re-sealable bag. Always KEEP the desiccant within the bag.

# CHAPTER 2

## SAMPLE PREPARATION





## 2.1. Sample Requirements for Fresh Frozen Tissue




To prevent RNA degradation, we recommend performing tissue embedding within **30 min** upon harvest.




This manual also supports frozen sample embedding. However, if not properly stored, frozen samples may have fragmentation, deformation, or other suboptimal conditions. In addition to RNA integrity, tissue morphology should be carefully assessed before Stereo-seq experiments. It is crucial to ensure that RNA integrity and tissue morphology both meet the transcriptomics requirements before starting.

Based on chip specifications, there are three tissue sizes:


1. Stereo-seq chip (0.5 cm \* 0.5 cm): The tissue size should not exceed **0.45 cm x 0.45 cm x 2 cm**, considering that the tissue section should not exceed 80% area coverage of the chip.



2. Stereo-seq chip (1 cm \* 1 cm): The tissue size should not exceed **0.9 cm x 0.9 cm x 2 cm**, considering that the tissue section should not exceed 80% area coverage of the chip.



3. Stereo-seq chip (1 cm \* 2 cm): The tissue size should not exceed **0.9 cm x 1.8 cm x 0.7 cm**;  
Stereo-seq chip (2 cm \* 2 cm): The tissue size should not exceed **1.8 cm x 1.8 cm x 0.7 cm**;  
Stereo-seq chip (2 cm \* 3 cm): The tissue size should not exceed **1.8 cm x 2.7 cm x 0.7 cm**;  
considering that the tissue section should not exceed 80% area coverage of the chip.



## Sample Types

For Stereo-seq Transcriptome FF Solution, Stereo-seq Transcriptomics mIF Solution, and Stereo-seq Large Chip Designs Solution: These product solutions are compatible with samples from all common animals, including but not limited to human, monkey, and mouse.

For Stereo-CITE Proteo-Transcriptomics Solution: This product solution is compatible with the simultaneous detection of multiple proteins and whole transcriptomes in both human and mouse samples.

For details, refer to the STOmics Validated Tissue List at:

<https://en.stomics.tech/resource/STOmicsTestedTissueList?lang=en#>

## FF Sample RNA Integrity Number (RIN) Value

It is recommended that you check the RNA quality (RIN value) of a tissue sample before proceeding to the Stereo-seq experiment. Collect 10 - 20 sections of 10  $\mu$ m and store them at -20°C in a pre-cooled 1.5 mL EP tube for total RNA extraction. Refer to Figure 1 for the peak RNA RIN value in mouse brain tissue sections.



1. If you are using Stereo-seq Transcriptomics FF V1.3 Solutions, Stereo-seq Transcriptomics FF V1.3 for mIF Solutions, and Stereo-CITE Proteo-Transcriptomics V1.1 Solutions, it is recommended that you perform subsequent experiment operations only on tissue samples with **RIN value  $\geq 4$** .
2. If you are using Stereo-seq Transcriptomics Large Chip Designs, it is recommended that you perform subsequent experiments only on tissue samples with **RIN value  $\geq 6$** .
3. If you are using Stereo-seq Transcriptomics FF V1.2 Solutions, Stereo-seq Transcriptomics FF V1.2 for mIF Solutions, Stereo-CITE Proteo-Transcriptomics V1.0 Solutions, and Stereo-seq Transcriptomics Large-size Chip, it is recommended that you conduct subsequent experiments solely on tissue samples with **RIN value  $\geq 7$** .




Figure 1. Example of RNA size distribution and RIN value measurement of tissue sections.

## 2.2. Sample Embedding for Fresh Frozen Tissue

Prepare the equipment and materials in advance. The following are the consumables or reagents required for embedding one sample. If multiple samples need to be embedded, use more consumables or reagents as necessary.



Materials		
Brand	Description	Quantity
-	Crushed ice in a container	1
-	Dry ice in a container	1
-	Aluminum Foil	1
-	Sealable Plastic Bag	1
BIOSHARP/Metal Coolbox/BC032	Metal Block	1
-	Sterile Gauze	2
Corning	Corning® 35 mm TC-treated Culture Dish (353001)	1
Sakura/Base Molds/4583	O.C.T. Compound	1
Sakura/Base Molds/4132/4133 (or other brands of appropriately-sized embedding molds)	Stainless-steel Base Mold A	1
Sakura/Base Molds/4133/41565/4124 (or other brands of appropriately-sized embedding molds)	Stainless-steel Base Mold B (slightly larger than mold A)	1

## Materials

Brand	Description	Quantity
-	Blunt-end Forceps	1
-	Syringe	1
-	Spatula	1
-	Scissors	1
-	Steel Ruler	1



For a demonstration video of tissue embedding, refer to the link below or scan the QR code:


<https://en.stomics.tech/resources/videos>

a. Prepare these materials in advance:


1. Pre-cool OCT on ice for **10 min** in advance.
2. Use **2** stainless-steel base molds slightly larger than the tissue of interest: **mold A** and **mold B** (slightly larger than **mold A**).
3. Add a few drops of pre-cooled OCT in **mold A** until it reaches approximately 2/3 of the mold and pre-cool on ice for **> 10 min** (remove introduced bubbles using a syringe).
4. In a petri dish filled with OCT, pre-cool the OCT on ice for **> 10 min** (remove introduced bubbles using a syringe).



5. For this step, you will need a box of dry ice and a metal block that has a flat surface to support the stainless-steel base mold when placed on dry ice. The size of the metal block should be larger than the stainless-steel base mold. Place the metal block on dry ice and pre-cool for **> 5 min** with the flat surface facing up.



6. Place **mold B** on dry ice and pre-cool for **> 5 min** .



b. Upon harvesting within **30 min**, use sterile gauze or dust-free paper to absorb excess liquid on the tissue surface to prevent ice formation in later steps.



c. Place the tissue in pre-cooled OCT and wrap the tissue evenly with OCT using a spatula without introducing air bubbles.



d. Remove any air bubbles using a syringe.







e. Orient the tissue such that the side intended for sectioning faces down, and then place it into **mold A**. Make sure the tissue is at the bottom of **mold A**, and then fill the mold with chilled OCT, without introducing bubbles, until the tissue is fully covered (remove any air bubbles using a syringe).



f. Place the tissue containing **mold A** onto the metal block that was placed on dry ice.



g. Place the pre-cooled steel ruler on the long and high side of **mold A** (to prevent the tissue from being compressed and deformed). Use **mold B** as a lid with the opening facing up, gently place it on top of **mold A**, and then place a few dry ice cubes on top of **mold B**. Make sure the two stainless-steel base molds can be completely covered with dry ice cubes.



h. After **5 min**, remove **mold B** and the steel ruler. Confirm that the OCT is completely frozen and has turned opaque; otherwise, repeat step f.



i. If the tissue block has solidified and turned opaque, grip the two edges of **mold A** and press the edges down to detach the tissue block from the mold.



j. Check the sectioning side of the tissue and determine if it is completely covered with OCT. If it is not completely covered, place the tissue block on the metal block, sectioning side up, add a few drops of OCT, and then wait until it solidifies and turns opaque.



k. Label the tissue block to mark the orientation of the tissue.



## 2.3. Sample Storage and Transportation

For storing, wrap the tissue block with aluminum foil and keep it in a properly labeled sealable plastic bag to prevent dehydration and damage, then store at  $-80^{\circ}\text{C}$ . For transportation, ship samples **on dry ice** in accordance with local regulations.