

Stereo-CITE Proteo-Transcriptomics Set

User Guidance

01 Introduction

STOmics Stereo-CITE Proteo-Transcriptomics Set is intended to simultaneously detect the whole transcriptome and high-plex protein on the same tissue section. Built upon DNA Nanoball (DNB) technology, STOmics Stereo-CITE Proteo-Transcriptomics Set enables a "tissue-to-data" solution through in situ capture of the whole transcriptome and high-plex protein markers at nanoscale resolution and centimeter-sized field of view (FOV). The Stereo-seq Chip T (poly-T-based chip) is loaded with capture probes containing spatial coordinate information. Through a series of biochemical processes, the probes can capture mRNA molecules and antibody-derived tags (ADTs) *in situ* within the tissue, then through cDNA synthesis, and obtain transcriptome plus multi-protein spatial distribution information of the entire tissue through sequencing and a complementary visualization platform.

02 Workflow



Figure 1 Stereo-CITE Proteo-Transcriptomics Set V1.1 Workflow





Figure 2 Overview of Stereo-CITE Proteo-Transcriptomics Set V1.1 Data Analysis Workflow

03 Resources Index

Step	Document Title and Description	Link
Experiment Preparation	Stereo-CITE Proteo-Transcriptomics Set User Guidance Document No. : STUM-UG005 This guidance aims to give the user a general view of the Stereo-CITE Proteo-Transcriptomics Set V1.1, with the overview of whole experiment workflow and index of resources to support the users' in-house experiments.	https://en.stomics.tech/re-sources/docu-ments/list.html 
	Microscope Assessment Guideline Document No. : STUM-PE001 This manual aims to guide the users to determine a proper microscope for STOmics application, introducing the microscope hardware requirements as well as guidance of imaging acquisition and evaluation.	
	Stereo-seq Chip Slide Operation Guide for Receiving, Handling and Storing This manual provides receiving, handling and storing guidance of the Stereo-seq Chip Slide.	
	Experiment Checklist For Stereo-CITE Proteo-Transcriptomics Set V1.1 This manual provides a detailed checklist for Stereo-CITE Proteo-Transcriptomics Set V1.1 workflow in terms of different experiment stages.	
	Sample Preparation Guide for Fresh Frozen Samples on Stereo-seq Chip Slides Document No. : STUM-SP001 This manual aims to guide the users in fresh frozen samples preparing for Stereo-seq Transcriptomics Solution V1.3 for mIF applications.	
Permeabilization Optimization and Transcriptomics Workflow	Stereo-seq Permeabilization Set for Stereo-CITE Proteo-Transcriptomics Application User Manual Document No. : STUM-PR004 This manual aims to guide the users in permeabilization optimization for specific samples that are used for Stereo-CITE Proteo-Transcriptomics applications.	
	Stereo-CITE Proteo-Transcriptomics Set User Manual Document No. : STUM-TT005 This manual provides a standard operation guidance for Stereo-CITE Proteo-Transcriptomics applications. This manual is suitable for Stereo-seq Chip T Slide(1cm * 1cm).	
Library Preparation and Sequencing	Stereo-CITE Proteo-transcriptomics Library Preparation User Manual Document No. : STUM-LP003 This manual aims to provide guidance for the whole-transcriptome library construction from cDNA and ADT products obtained via Stereo-CITE Proteo-transcriptomics workflow.	https://www.-completegenom-ics.com/docu-mentation/
	Related Sequencing Manuals CG DNBSEQ-T7RS Stereo-seq Visualization Reagent Set Instructions for Use CG DNBSEQ-G400RS Stereo-seq Visualization Reagent Set Instructions for Use	
	MGI DNBSEQ-T7RS Stereo-seq Visualization Reagent Set Instructions for Use MGI DNBSEQ-G400RS Stereo-seq Visualization Reagent Set Instructions for Use	
Data Analysis	StereoMap User Manual StereoMap is a desktop application designed to provide the essential analysis functionality you need to explore your Stereo-seq data interactively. This manual provides guidance in using StereoMap.	https://en.stomics.tech/service/-stereoMap-operation-manual.html
	SAW User Manual Stereo-seq Analysis Workflow (SAW) software suite is a set of pipelines bundled to map sequenced reads to their spatial location on the tissue section, quantify spatial feature expression, and visually present spatial expression distribution. SAW processes the data from the Stereo-seq sequencing platform, combined with microscope images, to generate spatial feature expression matrices. Analysts can use the output files as a starting point to perform downstream analysis. This manual provides guidance in using SAW.	https://en.stomics.tech/service/new-saw-operation-manual.html

04 Revision History

Version: A Date: Mar. 2025 Description: Initial release	<p>©2025 STOmics Tech Co., Ltd. All rights reserved.</p> <ol style="list-style-type: none">The products shall be for research use only, not for use in diagnostic procedures.The contents of this manual may be protected in whole or in part by applicable intellectual property laws. STOmics Tech Co., Ltd. and/or corresponding right subjects own their intellectual property rights according to law, including but not limited to trademark rights, copyrights, etc.STOmics Tech Co., Ltd. does not grant or imply the right or license to use any copyrighted content or trademark (registered or unregistered) of ours or any third party's. Without our written consent, no one shall use, modify, copy, publicly disseminate, change, distribute, or publish the program or contents of this manual without authorization, and shall not use the design or the design skills to use or take possession of the trademarks, the logo, or other proprietary information (including images, text, web design or form) of ours or those of our affiliates.Nothing contained herein is intended to or shall be construed as any warranty, expression or implication of the performance of any products listed or described herein. Any and all warranties applicable to any products listed herein are set forth in the applicable terms and conditions of sale accompanying the purchase of such product. STOmics Tech Co., Ltd. makes no warranty and hereby disclaims any and all warranties as to the use of any third-party products or protocols described herein.
---	--

05 Contact Us



STOmics Website: <https://en.stomics.tech>