

# **STOmics**

**Stereo-seq Chip P Slide (1 cm \* 1 cm)**

**Stereo-seq Chip T Slide (1 cm \* 1 cm)**

**Operation Guide For Receiving  
Handling And Storing**

## Upon Receiving:

- Stereo-seq Chip Slides are shipped at room temperature in vacuum sealed aluminum bags.
- Upon receiving, immediately store unopened Stereo-seq Chip Slides at  $-20^{\circ}\text{C}$  or  $4^{\circ}\text{C}$ .
- It can be store up to the expiration date on the label at  $-20^{\circ}\text{C}$  or  $4^{\circ}\text{C}$  unopened.

## Once Opened:

- Place unused Stereo-seq Slide back in the slide container and in the original aluminum bag. Keep sealed with the desiccant package inside the aluminum bag.
- Store at  $-20^{\circ}\text{C}$  or  $4^{\circ}\text{C}$  .